

World Wide Exclusive Formula

New **bluefusion**TM Gel Technology

Advanced Air
Pocket Technology,
32 air pockets per inch

The worlds first pre-loaded whitening
tray with blue fusion technology.
Designed by World Leading Dentist Dr G.
Favero & CEO Bryan J Christie

Patent pending

Blue Catalyst
Activator with
Photo Initiator

Ergonomically designed
with a greater depth for
maximum coverage

Soft comfortable and flexible
Megawhite tray, one-size fits all

Made from the highest quality
food gradeable materials

Self Service Instructions

1



You open the
bag at slits

2



You bite into
the gel tray

3



You position the lamp
in front of your mouth
then press play twice

4



You remove the mouth tray
after 15 minutes and place in
the hygiene bag provided